



DATE: November 20, 2015

TO: All physicians, Residents, Clinical Managers and Educators at SJH and HHSC

CC: Tracy Carrier, Albert Incretolli, Karon Taggart

FROM: Dr. P. Kavsak, Clinical Chemistry and Immunology

RE: **New Acetaminophen Assay in Hamilton**

Please be advised that effective December 1, 2015, we will be implementing a new Acetaminophen test (the Acetaminophen L3K® reagent from Sekisui Diagnostics) measured on the Abbott ARCHITECT platforms (i.e., all Hamilton Hospitals). Importantly, this assay is not subject to interference by NAC (N-acetylcysteine).

Accordingly, we will be removing the following comment on test results generated with this new reagent: *“Acetaminophen levels may be falsely low immediately after administration of N-acetylcysteine”*.

However, during our evaluation of NAC interference we have identified other tests that are affected by the presence of NAC¹.

Please note if a patient is receiving NAC, then lipase, lactate and possibly uric acid, total cholesterol and enzymatic creatinine (CRNEO – Creatinine Neonatal) may be falsely lower.

For further information please contact:

Dr. Tony Chetty
Discipline Director Clinical Chemistry and Immunology
905-522-1155, ext 34366
chetty@hhsc.ca

Dr. Peter Kavsak
Clinical Chemistry and Immunology
905-521-2100, ext 76293
kavsak@hhsc.ca

1. Agbor TA, Brown M, Clark L, Kavsak PA. Negative interference of N-acetyl cysteine (NAC) on selected chemistries on the Abbott architect platform. Clin Chim Acta. 2015 Oct 3.